



NITRILE EXAM GLOVES

POWDER-FREE



information and declaration
of conformity according to
PPER(EU) 2016/425
<http://www.kynex.nl/certificates>

FEATURES

- Excellent softness and wearing fitness
- Finger textured for a reliable, consistent grip
- Beaded cuff ensures easy donning and prevents roll down
- Compliant with food handling requirements
- Latex-free, powder-free and odor-free
- Sustainable choice: Due to the efficient packaging the CO₂ emission is reduced



PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Available Sizes	XS,S,M,L,XL,XXL	Color	Black
Glove Length	≥240mm	Freedom from Holes (Inspection level I)	AQL 1.5
Powder Content	Powder-Free	Sterile	No
Finger Thickness(mm)	0.10±0.03mm	Palm Thickness (mm)	0.08±0.03mm
Force at break	≥6N before and after aging	Tensile strength	≥14MPa before and after aging
Protection against Bacteria and Fungi	pass	Protection against 40% Sodium Hydroxide (K) 30% Hydrogen peroxide (P) 37% Formaldehyde (T)	> 30min
Protection against Viruses	pass		
Dimensions(±10mm)	XS/70;S/80;M/95;L/110;XL/120;XXL/120	Surface	Micro-textured on fingers

ORDERING INFORMATION

Color	Gram	Size						Type
		XS	S	M	L	XL	XXL	
Black	3.5	8403506	8403507	8403508	8403509	8403510	8403511	MDR+PPE

PACKAGING AND STORAGE

Packaging	100pcs/box,10boxes/ctn, 800 boxes/pallet, pallet dimension (80 x 120 x 195 cm)
Box size	21,5 x 12 x 6 cm
Carton Dimension	31,5 x 25 x 22,5 cm
Shelf Life	5 Years (recommended)
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

Compliant with EN455-1:2020, EN455-2:2015, EN455-3:2015, EN455-4:2009
Compliant with EN ISO 374-1:2016+A1:2018, EN ISO 374-5:2016
Notified body CE2777 (SATRA Technology Europe Ltd)